
## To which Arduino pins does the Vernier Interface Shield conne

MAY 12, 2017

The <u>Vernier Arduino</u>® <u>Interface Shield</u> (BT-ARD), has four sensor ports; two analog and two dig the following Arduino pins:



## **Analog 1**

Sensor pin 1 (+/-10 V output line) — Arduino pin A1

Sensor pin 2 (Ground) – Arduino ground

Sensor pin 3 (Vres/resistance reference) – Arduino pin A4

(also for I2C autoID)

Sensor pin 4 (AutoID) – Arduino pin A5 (also for I2C autoID)

Sensor pin 5 (+5 Volt power) – Ardino 5V

Sensor pin 6 (0-5 V output line) – Arduino pin A0

## Digital 2

Sensor pin 1 (Input/Output Line 0) – Arduino pin D6

Sensor pin 2 (Input/Output Line 1) – Arduino pin D7/A4

Sensor pin 3 (Input/Output Line 2) – Arduino pin D8/A5

Sensor pin 4 (+5 Volt power) – Ardino 5V

Sensor pin 5 (Ground) – Arduino ground

Sensor pin 6 (Input/Output Line 3) – Arduino pin D9

## Notes:

- \*Details for Vernier sensor pinouts can be found here: <a href="https://www.vernier.com/support/sensor-p">https://www.vernier.com/support/sensor-p</a>
- \*Most analog sensors use the 0-5 V output line.
- \*Digital sensors vary in how they use the four IO lines.
- \*Vernier interfaces use sensor pin 3 and 4 to auto-ID connected sensors. These pins are connected t code for auto-ID'ing sensors can be found <a href="here">here</a>. For many modern Vernier sensors, lines A4 and A sensors. On some Arduino types, notably Arduino Mega, these pins cannot be used for I2C commu Arduino Mega uses different pins for I2 communications.
- \*Pins A4 and A5 are shared across all four connectors. In order to properly access the BTA and BT the analog multiplexer circuit; MUX Control addresses are D10 and D11.
- \*An LED indicator is tied to D13 and a general purpose button is connected to D12.

For more information about how Vernier sensors work with Arduino, see our sample code and pro